

# plum QuickFix

# FIX



## Your quick help for minor everyday accidents

Mount the plaster dispenser on the wall and save time.  
Put on the plaster - and back to work!

### QUICK AND HYGIENIC HELP WITH MINOR INJURIES

Stop wasting time on minor injuries. Grab a plaster from the dispenser, put it on, and go back to work! Plum QuickFix solves the problem easily, hygienically and effectively.

Lots of minor accidents happen every day at work where a plaster is all that is needed. Choose the right plaster for your work situation. A suitable solution is available whether you need plasters for common work situations, for workplaces with a lot of humidity or for work within the food sector.



### PLUM QUICKFIX IS A UNIQUE DISPENSER

- Immediate and visible access to plasters and eyewash ampoules
- Each plaster is ready for use with one pull
- Apply with one hand
- Easily refilled as required
- Easy to see when it is empty
- Quickly stops minor bleeding
- Small dimensions - easy to fit

# Flexible system for all types of work

- Plum QuickFix provides a structured approach to first aid of minor injuries



## PLUM QUICKFIX DISPENSER

The classic green dispenser contains two refills of Plum QuickFix plasters and/or Plum Eyewash 0.9% Sodium Chloride QuickRinse ampoules. The green cap protects the content of the dispenser from dust and dirt from the workplace. It is designed with an inside shelf for smaller objects and it can be secured with lock and key.

The dispenser can be used in any location with risk of injuries and scratches. The dispenser is designed for mounting on any stationary workplace.

Dimensions: h: 135 x w: 230 x d: 32 mm



Plum QuickFix Uno, white



Plum QuickFix Uno, transparent



Plum QuickFix Solo®, transparent

## PLUM QUICKFIX UNO AND PLUM QUICKFIX SOLO® DISPENSERS

For workplaces that are not challenged with dust and dirt, we have developed the Plum QuickFix Uno and Plum QuickFix Solo® plaster dispensers. They are small and modern dispensers in a unique design that combines workplace safety with modern style.

Plum QuickFix Uno:  
Dimensions: h: 85 x w: 130 x d: 35 mm

Plum QuickFix Solo® :  
Dimensions: h: 80 x w: 130 x d: 27 mm



# How to use the plaster dispenser

## Step 1



Grab a plaster

## Step 2



Put it on with one hand

## Step 3



And back to work

# Choose your Plum QuickFix plaster dispenser

- to meet your specific needs

Lots of minor accidents happen every day at work where a plaster is all that is needed



**PLUM QUICKFIX WATER RESISTANT**  
Dispenser with 90 water resistant plasters suitable for environments with a high level of humidity.

Packaging: 10  
Article no: 5501  
Refill article no: 5511



**PLUM QUICKFIX BLOOD STOPPER**  
Dispenser with 90 blood stopper plasters with hemostatic effect.

Packaging: 10  
Article no: 5510  
Refill article no: 5516



**PLUM QUICKFIX&RINSE**  
Dispenser with 45 elastic plasters and 5 ampoules with 20 ml 0.9% sterile sodium chloride.

Packaging: 10  
Article no: 4631  
Refill article no: 5512 & 5160



**PLUM QUICKFIX BLACK**  
Dispenser with 90 black plasters made from a durable, black polyester fabric.

Packaging: 10  
Article no: 5536  
Refill article no: 5522



**PLUM QUICKFIX ELASTIC**  
Dispenser with 90 elastic plasters that adjust to the movement of the skin.

Packaging: 10  
Article no: 5502  
Refill article no: 5512



**PLUM QUICKFIX ELASTIC & WATER RESISTANT**  
Dispenser with 45 elastic plasters and 45 water resistant plasters.

Packaging: 10  
Article no: 5507  
Refill article no: 5512 & 5511



**PLUM QUICKRINSE**  
Dispenser with 10 eyewash ampoules with 20 ml 0.9% sterile sodium chloride.

Packaging: 10  
Article no: 4632  
Refill article no: 5160



**PLUM QUICKFIX BLACK LONG**  
Dispenser with 60 long black plasters made from a durable, black polyester fabric.

Packaging: 10  
Article no: 5538  
Refill article no: 5523



**PLUM QUICKFIX DETECTABLE**  
Dispenser with 90 blue plasters with built-in metal surface which can be traced by a metal detector.

Packaging: 10  
Article no: 5503  
Refill article no: 5513



**PLUM QUICKFIX ELASTIC LONG**  
Dispenser with 60 long elastic plasters that can be wrapped around the finger several times.

Packaging: 10  
Article no: 5528  
Refill article no: 5508



**PLUM QUICKFIX XRAY**  
Dispenser with 90 XRay plasters that are water resistant and metal and X-ray detectable.

Packaging: 10  
Article no: 5537  
Refill article no: 5524



**PLUM QUICKFIX XRAY LONG**  
Dispenser with 60 long XRay plaster that are water resistant and X-ray detectable.

Packaging: 10  
Article no: 5539  
Refill article no: 5525



**PLUM QUICKFIX ALOE VERA**  
Dispenser with 90 aloe vera plasters, which soothes and moisturizes.

Packaging: 10  
Article no: 5506  
Refill article no: 5514



**PLUM QUICKFIX DETECTABLE LONG**  
Dispenser with 60 long blue plasters with built-in metal surface which can be traced by a metal detector.

Packaging: 10  
Article no: 5529  
Refill article no: 5509



**PLUM QUICKFIX EMPTY DISPENSER**  
Dispenser without refills. Fill with any Plum QuickFix or Plum QuickRinse refills of your choice.

Packaging: 10  
Article no: 5500

# Plum QuickFix Uno and Solo®

Fill the dispenser with the plaster of your own choosing. All the QuickFix plaster refills can be used according to your specific needs.

The dispenser can also be used with a Plum QuickRinse refill.



**PLUM QUICKFIX UNO, TRANSPARENT**  
Plum QuickFix Uno, transparent dispenser without refills.

Article no: 5533  
Packaging: 10



**PLUM QUICKFIX UNO ELASTIC, TRANSPARENT**  
Transparent dispenser with 45 elastic and comfortable plasters that adjust to the movement of the skin. The plasters are skin-friendly and allow the skin to breathe.

Article no: 5531  
Packaging: 10



**PLUM QUICKFIX UNO, WHITE**  
Plum QuickFix Uno, white dispenser without refills.

Article no: 5534  
Packaging: 10



**PLUM QUICKFIX UNO ELASTIC, WHITE**  
White dispenser with 45 elastic and comfortable plasters that adjust to the movement of the skin. The plasters are skin-friendly and allow the skin to breathe.

Article no: 5532  
Packaging: 10



PLUM QUICKFIX SOLO® Transparent - shown with example of plaster refill



**PLUM QUICKFIX SOLO®, TRANSPARENT**  
Plum QuickFix Solo®, transparent dispenser without refills.

Article no: 5535  
Packaging: 2

# How to fill the Plum QuickFix dispenser

## Step 1



Remove the bottom part of the refill cardboard

## Step 2



Push the refill pack into the dispenser until you hear a 'click'

## Step 3



The plasters are now ready for use

# Plum QuickFix plaster refills



**PLUM QUICKFIX ELASTIC**  
45 elastic and comfortable textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe.

**REFILL**  
Display carton: 6 x 45  
Packaging: 8 x 6 x 45  
Article no: 5512  
Dimensions: l: 72 x w: 25 mm

**DISPENSER**  
Article no: 5502



**PLUM QUICKFIX DETECTABLE**  
45 blue conspicuous textile plasters with a built-in metal foil that can be traced by a metal detector. Particularly suitable for the food industry. The plasters are made from a breathable material.

**REFILL**  
Display carton: 6 x 45  
Packaging: 8 x 6 x 45  
Article no: 5513  
Dimensions: l: 72 x w: 25 mm

**DISPENSER**  
Article no: 5503



**PLUM QUICKFIX ALOE VERA**  
The aloe vera contained in the pad has well known soothing, moisturizing and calming properties. It creates the best conditions for the wound environment helping the natural healing process. Prevents the formation of scabs. The material of the plaster is textile and breathable.

**REFILL**  
Display carton: 6 x 45  
Packaging: 8 x 6 x 45  
Article no: 5514  
Dimensions: l: 72 x w: 25 mm

**DISPENSER**  
Article no: 5506



**PLUM QUICKFIX WATER RESISTANT**  
45 water resistant plasters that are particularly suitable for working environments with a high level of humidity. The plasters are made from a protective PE material and is breathable.

**REFILL**  
Display carton: 6 x 45  
Packaging: 8 x 6 x 45  
Article no: 5511  
Dimensions: l: 72 x w: 25 mm

**DISPENSER**  
Article no: 5501



**PLUM QUICKFIX BLOOD STOPPER**  
Plasters with hemostatic effect. The presence of alginate fibres between the nonwoven fibres of the wound pad creates an artificial clot (gel) inside the pad, which blocks the bleeding in 1-2 minutes and maintains optimal moist conditions. The plaster is breathable.

**REFILL**  
Display carton: 6 x 45  
Packaging: 8 x 6 x 45  
Article no: 5516  
Dimensions: l: 72 x w: 25 mm

**DISPENSER**  
Article no: 5510



**PLUM QUICKFIX BLACK**  
This pack contains 45 plasters made from black polyester fabric. The plasters are designed to be durable and breathable, making them suitable for various workplaces.

**REFILL**  
Display carton: 6 x 45  
Packaging: 8 x 6 x 45  
Article no: 5522  
Dimensions: l: 72 x w: 25 mm

**DISPENSER**  
Article no: 5536



**PLUM QUICKFIX XRAY**  
45 blue textile X-ray detectable plasters with an aluminum inlay containing barium sulphate. The plasters are water resistant and detectable by metal detectors and X-rays.

**REFILL**  
Display carton: 6 x 45  
Packaging: 8 x 6 x 45  
Article no: 5524  
Dimensions: l: 72 x w: 25 mm

**DISPENSER**  
Article no: 5537

	Elastic (*1)	Detectable (*1)	Water Resistant (*1)	Aloe Vera (*2)	Blood Stopper (*3)	Black (*1)	XRay (*1)
<b>Material</b>	Pink cross direction elastic fabric (Cotton & polyamide)	Blue cross direction elastic fabric (Cotton & polyamide). Detectable aluminium foil facing the backing.	Pink perforated PE (Polyethylene)	Pink cross direction elastic fabric (65% cotton & 35% polyamide)	Pink cross direction elastic fabric (65% cotton & 35% polyamide)	Black Polyester fabric	Blue Ethylene Vinyl Acetate and PE (Polyethylene) fabric with an aluminium foil inlay containing barium sulphate
<b>Adhesive pad</b>	Hypoallergenic acrylic polymer	Hypoallergenic acrylic polymer	Hypoallergenic acrylic polymer	Hypoallergenic acrylic polymer	Hypoallergenic acrylic polymer	Hypoallergenic hot melt polymer	Hypoallergenic acrylic polymer
<b>Absorbing pad</b>	Viscose and non-woven PP	Viscose and PP/PE non-woven	Viscose and non-woven PP	Viscose and PP non-woven with 1,5% aloe vera	Viscose and PP/PE non-woven with haemostatic layer (calcium-sodium alginate particles) and anti-adhering layer (polyester net)	Viscose and polyester	Viscose and polyester
<b>Special Feature</b>	- Elastic - Breathable	- Elastic - Metal detectable - Breathable	- Water resistant - Breathable	- Soothing and moisturizing - Breathable	- Haemostatic effect - Breathable	- Breathable	- Water resistant - X-ray and metal detectable

\*1: According to TDS C1 version 2.0, 08-08-2025

\*2: According to TDS C2a Aloe Vera version 2.0, 08-08-2025

\*3: According to TDS C2a Blood Stopper version 2.0, 08-08-2025

# Plum QuickFix plaster refills - long



**PLUM QUICKFIX ELASTIC LONG**  
 Display carton: 6 x 30  
 Packaging: 8 x 6 x 30  
 Article no: 5508  
 Dimensions: l: 120 x w: 20 mm

**DISPENSER**  
 Article no: 5528

**PLUM DETECTABLE LONG**  
 Display carton: 6 x 30  
 Packaging: 8 x 6 x 30  
 Article no: 5509  
 Dimensions: l: 120 x w: 20 mm

**DISPENSER**  
 Article no: 5529

**PLUM QUICKFIX LONG**  
 30 textile plasters with extra length that can be wrapped around the finger several times. Choose between elastic plasters, blue detectable plasters, black plasters, or x-ray plasters.



**PLUM QUICKFIX BLACK LONG**  
 Display carton: 6 x 30  
 Packaging: 8 x 6 x 30  
 Article no: 5523  
 Dimensions: l: 120 x w: 20 mm

**DISPENSER**  
 Article no: 5538

**PLUM XRAY LONG**  
 Display carton: 6 x 30  
 Packaging: 8 x 6 x 30  
 Article no: 5525  
 Dimensions: l: 120 x w: 20 mm

**DISPENSER**  
 Article no: 5539

# Plum QuickRinse

- for smaller accidents

Plum Eyewash 0.9% Sodium Chloride QuickRinse is a sterile first aid eye irrigation.

The product can be used for washing away:

- Dust
- Dirt
- Foreign bodies

Plum Eyewash 0.9% Sodium Chloride QuickRinse ampoules are designed specifically for rinsing away small foreign particles that can be removed with just a small amount of fluid. A key advantage of these ampoules is their convenience – you only need to open a small, single-use container rather than a larger bottle.




**PLUM EYEWASH 0.9% SODIUM CHLORIDE QUICKRINSE - 5 X 20 ML**  
 5 eyewash ampoules with 20 ml 0.9% sterile sodium chloride. For use individually or for refilling Plum QuickSafe®, Plum QuickFix&Rinse and Plum QuickRinse dispensers.

**Content:** 5 x 20 ml  
**Packaging:** 4 x 5 x 20 ml  
**Article no:** 5160

## INSTRUCTIONS FOR USE

Plum Eyewash 0.9% Sodium Chloride QuickRinse is easy to use. Follow these simple steps:



**Step 1:**  
 Twist the ampoule



**Step 2:**  
 Rinse the eye



**Step 3:**  
 Throw away the empty ampoule



SAFETY AT WORK



Plum Safety ApS  
Mandelalleen 1  
DK-5610 Assens

T: +45 6916 9600  
info@plum.eu

[www.plum.eu](http://www.plum.eu)



info@plum.eu  
+45 6916 9600

[www.plum.eu](http://www.plum.eu)

**plum**  
Safety